


Product Datasheet

Human MICA ELISA Kit (orb390959)

Description	Sandwich ELISA for quantitative detection of human MICA concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	MICA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Cell Lysates, Tissue Homogenates, Serum
Sample Types	Cell Culture Supernates, Serum, Cell Lysates, Tissue Homogenates
Uniprot ID	Q29983
Sensitivity	10 pg/ml
Entrez	100507436
Expiration Date	Please enquire.



Representative
Standard
Curve