


Product Datasheet

Human MIF ELISA Kit (orb390963)

Description	Sandwich ELISA for quantitative detection of human MIF concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	MIF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P14174
Sensitivity	15 pg/ml
Entrez	4282
Expiration Date	Please enquire.



Representative
Standard
Curve