


Product Datasheet

Human MYD88 ELISA Kit (orb390975)

Description	Sandwich ELISA for quantitative detection of human MYD88 concentrations.
Reactivity	Human
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	MYD88
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants
Sample Types	Cell Culture Supernates
Uniprot ID	Q99836
Sensitivity	10 pg/ml
Entrez	4615
Expiration Date	Please enquire.



Representative
Standard
Curve