


Product Datasheet

Human NOV ELISA Kit (orb390981)

Description	Sandwich ELISA for quantitative detection of human NOV concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	NOV
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P48745
Sensitivity	2 pg/ml
Entrez	4856
Expiration Date	Please enquire.



Representative
Standard
Curve