


Product Datasheet

Human OPG ELISA Kit (orb390982)

Description	Sandwich ELISA for quantitative detection of human OPG concentrations.
Reactivity	Human
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	TNFRSF11B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	000300
Sensitivity	5 pg/ml
Entrez	4982
Expiration Date	Please enquire.



Representative
Standard
Curve