

Product Datasheet

Human P-Cadherin ELISA Kit (orb390983)

Description Sandwich ELISA for quantitative detection of human P-Cadherin concentrations.

Reactivity Human

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target CDH3

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum


Sample Types Cell Culture Supernates, Serum

Uniprot ID **P22223**

Sensitivity 2 pg/ml

Entrez **1001**

Expiration Date Please enquire.



Representative
Standard
Curve