


Product Datasheet

Human SPARC ELISA Kit (orb391000)

Description	Sandwich ELISA for quantitative detection of human SPARC concentrations.
Reactivity	Human
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	SPARC
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	P09486
Sensitivity	20 pg/ml
Entrez	6678
Expiration Date	Please enquire.



Representative
Standard
Curve