


Product Datasheet

Human TIMP1 ELISA Kit (orb391014)

Description	Sandwich ELISA for quantitative detection of human TIMP1 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	TIMP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P01033
Sensitivity	30 pg/ml
Entrez	7076
Expiration Date	Please enquire.



Representative
Standard
Curve