


Product Datasheet

Human FLT1 ELISA Kit (orb391026)

Description	Sandwich ELISA for quantitative detection of human FLT1 concentrations.
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	FLT1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P17948
Sensitivity	30 pg/ml
Entrez	2321
Expiration Date	Please enquire.



Representative
Standard
Curve