


Product Datasheet

Mouse BMP2 ELISA Kit (orb391035)

Description	Sandwich ELISA for quantitative detection of mouse BMP2 concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	BMP2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Bone Tissue Homogenates, Cell Culture Supernatants, Serum
Sample Types	Bone Tissue Homogenates, Cell Culture Supernates, Serum
Uniprot ID	P21274
Sensitivity	15 pg/ml
Entrez	12156
Expiration Date	Please enquire.



Representative
Standard
Curve