

Product Datasheet

Mouse TNFRSF8 ELISA Kit (orb391052)

Description Sandwich ELISA for quantitative detection of mouse TNFRSF8 concentrations.

Reactivity Mouse

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target TNFRSF8

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q60846**

Sensitivity 10 pg/ml

Entrez **21941**

Expiration Date Please enquire.



Representative
Standard
Curve