


Product Datasheet

Mouse CD40 ELISA Kit (orb391055)

Description	Sandwich ELISA for quantitative detection of mouse CD40 concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CD40
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Cell Lysates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	P27512
Sensitivity	1 pg/ml
Entrez	21939
Expiration Date	Please enquire.



Representative
Standard
Curve