


Product Datasheet

Mouse CLU ELISA Kit (orb391060)

Description	Sandwich ELISA for quantitative detection of mouse CLU concentrations.
Reactivity	Mouse
Range	18.7 ng/ml - 1200 ng/ml
Tested Applications	ELISA
Target	CLU
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q06890
Sensitivity	6 ng/ml
Entrez	12759
Expiration Date	Please enquire.



Representative
Standard
Curve