


Product Datasheet

Mouse DMP1 ELISA Kit (orb391071)

Description	Sandwich ELISA for quantitative detection of mouse DMP1 concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	DMP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	055188
Sensitivity	10 pg/ml
Entrez	13406
Expiration Date	Please enquire.



Representative
Standard
Curve