


Product Datasheet

Mouse FLT3LG ELISA Kit (orb391077)

Description	Sandwich ELISA for quantitative detection of mouse FLT3LG concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	FLT3LG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	P49772
Sensitivity	10 pg/ml
Entrez	14256
Expiration Date	Please enquire.



Representative
Standard
Curve