

Product Datasheet

IL18 antibody (Biotin) (orb39109)

Description	IL18 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Dog Interleukin-18 protein (37-193AA)
Target	IL18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9XSR0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.