

Product Datasheet

Mouse IL1 alpha ELISA Kit (orb391091)

Description Sandwich ELISA for quantitative detection of mouse IL1 alpha concentrations.

Reactivity Mouse

Range 1.9 pg/ml - 125 pg/ml

Tested Applications ELISA

Target IL1A

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID **P01582**

Sensitivity 1 pg/ml

Entrez **16175**

Expiration Date Please enquire.



Representative
Standard
Curve