

## Product Datasheet

### Mouse IL1 alpha ELISA Kit (orb391092)

**Description** Sandwich ELISA for quantitative detection of mouse IL1 alpha concentrations.

**Reactivity** Mouse

**Range** 1.9 pg/ml - 125 pg/ml

**Tested Applications** ELISA

**Target** IL1A

**Storage** Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


**Sample Types** Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

**Uniprot ID** **P01582**

**Sensitivity** 1 pg/ml

**Entrez** **16175**

**Expiration Date** Please enquire.



Representative  
Standard  
Curve