


Product Datasheet

Mouse IL21 ELISA Kit (orb391097)

Description	Sandwich ELISA for quantitative detection of mouse IL21 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	IL21
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	Q9ES17
Sensitivity	10 pg/ml
Entrez	60505
Expiration Date	Please enquire.



Representative
Standard
Curve