


Product Datasheet

Mouse LAMA1 ELISA Kit (orb391102)

Description	Sandwich ELISA for quantitative detection of mouse LAMA1 concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	LAMA1; LAMB1; LAMC1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P02469 , P02468 , P19137
Sensitivity	10 pg/ml
Entrez	16777 , 16772 , 226519
Expiration Date	Please enquire.



Representative
Standard
Curve