


Product Datasheet

Mouse MFGE8 ELISA Kit (orb391106)

Description	Sandwich ELISA for quantitative detection of mouse MFGE8 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	MFGE8
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P21956
Sensitivity	10 pg/ml
Entrez	17304
Expiration Date	Please enquire.



Representative
Standard
Curve