

Product Datasheet

Mouse Prolactin ELISA Kit (orb391119)

Description Sandwich ELISA for quantitative detection of mouse Prolactin concentrations.

Reactivity Mouse

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target PRL

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P06879**

Sensitivity 10 pg/ml

Entrez **19109**

Expiration Date Please enquire.



Representative
Standard
Curve