


Product Datasheet

Mouse SAA1 ELISA Kit (orb391125)

Description	Sandwich ELISA for quantitative detection of mouse SAA1 concentrations.
Reactivity	Mouse
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	SAA1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P05366
Sensitivity	800 pg/ml
Entrez	20208
Expiration Date	Please enquire.



Representative
Standard
Curve