


Product Datasheet

Mouse APCS ELISA Kit (orb391126)

Description	Sandwich ELISA for quantitative detection of mouse APCS concentrations.
Reactivity	Mouse
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	APCS
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P12246
Sensitivity	50 pg/ml
Entrez	20219
Expiration Date	Please enquire.



Representative
Standard
Curve