


Product Datasheet

Mouse TNF ELISA Kit (orb391133)

Description	Sandwich ELISA for quantitative detection of mouse TNF concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	TNF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P06804
Sensitivity	8 pg/ml
Entrez	21926
Expiration Date	Please enquire.



Representative
Standard
Curve