


Product Datasheet

Mouse TrkB ELISA Kit (orb391135)

Description	Sandwich ELISA for quantitative detection of mouse TrkB concentrations.
Reactivity	Mouse
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	NTRK2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Cell Lysates, Tissue Homogenates
Sample Types	Cell Culture Supernates, Cell Lysates, Tissue Homogenates
Uniprot ID	P15209
Sensitivity	5 pg/ml
Entrez	18212
Expiration Date	Please enquire.



Representative
Standard
Curve