


Product Datasheet

Mouse KDR ELISA Kit (orb391139)

Description	Sandwich ELISA for quantitative detection of mouse KDR concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	KDR
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P35918
Sensitivity	2 pg/ml
Entrez	16542
Expiration Date	Please enquire.



Representative
Standard
Curve