

Product Datasheet

PRL3 antibody (orb39114)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	PTP4A3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19535
Uniprot ID	O75365
NCBI	NP_009010.2 , NP_116000.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

