


Product Datasheet

Rat Cystatin C ELISA Kit (orb391149)

Description	Sandwich ELISA for quantitative detection of rat Cystatin C concentrations.
Reactivity	Rat
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	CST3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P14841
Sensitivity	10 pg/ml
Entrez	25307
Expiration Date	Please enquire.



Representative
Standard
Curve