

Product Datasheet

Rat IFN gamma ELISA Kit (orb391156)

Description Sandwich ELISA for quantitative detection of rat IFN gamma concentrations.

Reactivity Rat

Range 31.2 pg/ml - 2000 pg/ml

Tested Applications ELISA

Target IFNG

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates


Sample Types Cell Culture Supernates, Serum, Plasma, Tissue Homogenates

Uniprot ID [P01581](#)

Sensitivity 7 pg/ml

Entrez [25712](#)

Expiration Date Please enquire.



Representative
Standard
Curve