


## Product Datasheet

### Rat IGF2 ELISA Kit (orb391159)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of rat IGF2 concentrations.                         |
| <b>Reactivity</b>          | Rat                                                                                           |
| <b>Range</b>               | 46.9 pg/ml - 3000 pg/ml                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                         |
| <b>Target</b>              | IGF2                                                                                          |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                         |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Cell Culture Supernatants, Serum, Plasma      |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma                                                        |
| <b>Uniprot ID</b>          | <b>P01346</b>                                                                                 |
| <b>Sensitivity</b>         | 25 pg/ml                                                                                      |
| <b>Entrez</b>              | <b>24483</b>                                                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                               |



Representative  
Standard  
Curve