


Product Datasheet

Rat HAVCR1 ELISA Kit (orb391163)

Description	Sandwich ELISA for quantitative detection of rat HAVCR1 concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	HAVCR1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	054947
Sensitivity	2 pg/ml
Entrez	286934
Expiration Date	Please enquire.



Representative
Standard
Curve