


Product Datasheet

Rat NTF3 ELISA Kit (orb391167)

Description	Sandwich ELISA for quantitative detection of rat NTF3 concentrations.
Reactivity	Rat
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	NTF3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	P18280
Sensitivity	2 pg/ml
Entrez	81737
Expiration Date	Please enquire.



Representative
Standard
Curve