


Product Datasheet

Rat TGFB3 ELISA Kit (orb391173)

Description	Sandwich ELISA for quantitative detection of rat TGFB3 concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	TGFB3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	Q07258
Sensitivity	10 pg/ml
Entrez	25717
Expiration Date	Please enquire.



Representative
Standard
Curve