



Product Datasheet

MKP3 antibody (orb39120)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant Protein
Target	Dusp6
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42319
Uniprot ID	Q64346
NCBI	NP_446335.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...