

Product Datasheet

Human ADPRH protein (orb391205)

Description

Recombinant human ADPRH protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

ADPRH

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 50% Glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, 50% Glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SGLVPRGSH MEKYVAAMVL
SAAGDALGYY NGKWEFLQDG EKIHRQLAQL
GGLDALDVGR WRVSDDTVMH LATAEALVEA
GKAPKLTQLY YLLAKHYQDC MEDMDGRAPG
GASVHNAMQL KPGKPNGWRI PFNSHEGGCG
AAMRAMCIGL RFPHHSQLDL LIQVSIESGR
MTHHHPTGYL GALASALFTA YAVNSRPPLQ
WGKGLMELLP EAKKYIVQSG YFVEENLQHW
SYFQTKWENY LKLRGILDGE SAPTFPESFG
VKERDQFYTS LSYSWGWSGS GHDAPMIAYD
AVLAAGDSWK ELAHRFFHFG GDSDSTAAIA
GCWWGVMYGF KGVSPSNYEK LEYRNRLEET
ARALYSLGSK EDTVISL

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID
P54922