

Product Datasheet

Human APOH protein (orb391252)

Description

Recombinant human APOH protein is

system and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

APOH

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

GRTCPKPDDL PFSTVVPLKT FYEPGEEITY
SCKPGYVSRG GMRKFICPLT GLWPINTLKC
TPRVCPFAGI LENGAVRYTT FEYPNTISFS
CNTGFYLNKA DSAKCTEEGK WSPPLPVCAP
IICPPPSIPT FATLRVYKPS AGNNSLYRDT
AVFECLPQHA MFGNDTITCT THGNWTKLPE
CREVKCPFPS RPDNGFVNYP AKPTLYYKDK
ATFGCHDGYS LDGPPEIECT KLGNWSAMPS
CKASCKVPVK KATVVYQGER VKIQEKFKNG
MLHGDKVSFF CKNKEKKCSY TEDAQCIDGT
IEVPKCFKEH SSLAFWKTD A SDVKPCVDHH
HHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells