

Product Datasheet

Human BCL2L2 protein (orb391288)

Description

Recombinant human BCL2L2 protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

BCL2L2

Preservatives

Liquid in a 0.2 μM filtered solution of 25mM HEPES, 100mM KCl, 20% Glycerol, pH 7.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 25mM HEPES, 100mM KCl, 20% Glycerol, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MATPASAPDT RALVADFDVGY KLRQKGYVCG
AGPGEGPAAD PLHQAMRAAG DEFETRFRRRT
FSDLAAQLHV TPGSAQQRFT QVSDelfQGG
PNWGRLVAFF VFGAALCAES VNKEMEPLVG
QVQEWVMVAYL ETRLADWIHS SGGWAEFTAL
YGDGALEEAR RLREGNWASV RTLEHHHHHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q92843](#)

Entrez

[599](#)

Expiration Date

6 months from date of receipt.