

Product Datasheet

Human BLVRA protein (orb391300)

Description

Recombinant human BLVRA protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested
Applications

ELISA, MS, SDS-PAGE, WB

Target

BLVRA

Preservatives

Lyophilized from a 0.2 μM filtered solution of 4mM HCl.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 4mM HCl.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MERKFGVVVV GVGRAGSVRM RDLRNPHPSS
AFLNLIGFVS RRELGSIDGV QQISLEDALS
SQEVEVAYIC SESSSHEDYI RQFLNAGKHV
LVEYPM TSL AAAQELWELA EQKGKVSHEE
HVELLMEEFA FLKKEVVGKD LLKGSLLFTA
GPLEEERFGF PAFSGISRLT WLVS LFGELS
LVSATLEERK EDQYMKMTVC LETEKKSPLS
WIEEKGPGLK RNRYLSFHFK SGLENVPNV
GVNKNIFLKD QNIFVQKLLG QFSEKELAAE
KKRILHCLGL AEEIQKYCCS LEHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

P53004

Entrez

644