

Product Datasheet

Human C10ORF54 protein (orb391319)

Description

Recombinant human C10ORF54 protein is

system and the target gen...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

C10ORF54

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

FKVATPYSLY VCPEGQNVTL TCRLGPDVK
GHDVTFYKTW YRSSRGEVQT CSERRPIRNL
TFQDLHLHHG GHQAANTSHD LAQRHGLESA
SDHHGNFSIT MRNLTLDSG LYCCLVVEIR
HHHSEHRVHG AMELQVQTGK DAPSNVYYP
SSSQESENIT AAVDDIEGRM DEPKSCDKTH
TCPPCPAPEL LGGPSVFLFP PKPKDTLMIS
RTPEVTCVVV DVSHEDPEVK FNWYVDGVEV
HNAKTKPRE QYNSTYRVVS VLTVLHQDWL
NGKEYKCKVS NKALPAPIEK TISKAKGQPR
EPQVYTLPPS REEMTKNQVS LTCLVKGFYP
SDIAVEWESN GQPENNYKTT PPVLDSDGSF
FLYSKLTVDK SRWQQGNVFS CSVMHEALHN
HYTQKSLSLS PGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.