

Product Datasheet

Human C19ORF80 protein (orb391321)

Description

Recombinant human C19ORF80 protein is

system and the target gen...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

C19ORF80

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

CDKTHTCPPC PAPELLGGPS VLFPPKPKD
TLMISRTPEV TCVVVDVSHE DPEVKFNWYV
DGVEVHNAKT KPREEQYNST YRVVSVLTVL
HQDWLNGKEY KCKVSNKALP APIEKTISKA
KGQPREPQVY TLPPSRDELT KNQVSLTCLV
KGFYPSDIAV EWESNGQPEN NYKTTTPVLD
SDGSFFLYSK LTVDKSRWQQ GNVFSCSVMH
EALHNHYTQK SLSLSPGKDD DDKAPMGGPE
LAQHEELTLL FHGTLQLGQA LNGVYRTTEG
RLTKARNSLG LYGRTEILLG QEVSRGRDAA
QELRASLLET QMEEDILQLQ AEATAEVLGE
VAQAQKVL RD SVQRLEVQLR SAWLGPAYRE
FEVLKAHADK QSHILWALTG HVQRQRREMV
AQQHRLRQIQ ERLHTAALPA

Purity

Greater than 95% as determined by reducing SDS-PAGE.