

Product Datasheet

Mouse CADM1 protein (orb391341)

Description

Recombinant mouse CADM1 protein is

system and the target gene e...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

CADM1

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QNLFTKDVTV IEGEVATISC QVNKSDDSVI
QLLNPNRQTI YFRDFRPLKD SRFQLLNFS
SEKVS LTNV SISDEGRYFC QLYTDPPQES
YTTITVLVPP RNLMDIQKD TAVEGEEIEV
NCTAMASKPA TTIRWFKGNK ELKGKSEVEE
WSDMYTVTSQ LMLKVHKEDD GVPVICQVEH
PAVTGNLQTQ RYLEVQYKPKQ VHIQMTYPLQ
GLTREGDAFE LTCEAIGKPQ PVMVTWVRVD
DEMPQHAVLS GPNLFINN LN KTDNGTYRCE
ASNIVGKAHS DYMLYVYDPP TTIPPPTTTT
TTTTTTTTTI LTIITDTTAT TEPVHDSRA
GEEGTIGAVD HVDDIEGRMD EPKSCDKTHT
CPPCPAPELL GGPSVFLFPP KPKDTLMISR
TPEVTCVVVD VSHEDPEVKF NWYVDGVEVH
NAKTKPREEQ YNSTYRVVSV LTVLHQDWLN
GKEYKCKVSN KALPAIEKT ISKAKGQPRE
PQVYTLPPSR EEMTKNQVSL TCLVKGFYPS
DIKAEVWESNG OPENNYKTTT PVI DSDGSEF

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Source

Human cells