

Product Datasheet

Human CD99L2 protein (orb391437)

Description

Recombinant human CD99L2 protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

CD99L2

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DFDDFNLEDA VKETSSVKQP WDHTTTTTTN
RPGTTTRAPAK PPGSGLDLAD ALDDQDDGRR
KPGIGGRERW NHVTTTTTKRP VTTRAPANTL
GNDFDLADAL DDRNDRDDGR RKPIAGGGGF
SDKDLEDIVG GGEYKPDKGK GDGRYGSNDD
PGSGMVAEPG TIAVDDIEGR MDEPKSCDKT
HTCPPCPAPE LLGGPSVFLF PPKPKDTLMI
SRTPEVTCVV VDVSHEDPEV KFNWYVDGVE
VHNAKTKPRE EQYNSTYRVV SVLTVLHQDW
LNGKEYKCKV SNKALPAPIE KTISKAKGQP
REPQVYTLPP SREEMTKNQV SLTCLVKGFY
PSDIAVEWES NGQPENNYKT TPPVLDSDGS
FFLYSKLTVD KSRWQQGNVF SCSVMHEALH
NHYTQKSLSL SPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.