

Product Datasheet

Mouse CES2E protein (orb391469)

Description

Recombinant mouse CES2E protein is

system and the target gene e...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

CES2E

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MPLYKLLGWL NAVACGVLLL VLHVQQQDSA
SPIRNTHTGQ VRGSLVHVKD TDIHVHTFLG
IPFAKPPVGP LRFAPPEAPE PWSGVRDGT
HPNMCLQNDN LMGSEDLKMM NLILPPISMS
EDCLYLNIVY PAHAHEGSNL PVMVWIHGGA
LTVGMASMYD GSMLAATEDV VVVAIQYRLG
VLGFFSTGDQ HAKGNWGYLD QVAALRWVQQ
NIVHFGGNPD RVTIFGESAG GTSVSSHVVS
PMSQGLFHGA IMESGVAVLP DLISSSSEMV
HRIVANLSGC AAVNSETLMC CLRGKNEAEM
LAINKVFKII PGVVDGEFLP KHPQELMASK
DFHPVPSIIG INNDEYGWIL PTIMDPAQKI
EEITRKTLPV VLKSTALKMM LPPECGDLLM
EEYMGDTEDP ETLQAQFREM KGDFMFVIPA
LQVAHFQRSH APVYFYEFGH RPSFFKDFRP
PYVKADHGDE IFLVFGYQFG NIKLPYTEEE
EQLSRRIMKY WANFARHGPN NSEGLPYWVP
MDHDEFOYLOI DIOPSVGRAI KARRIOEFTK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558