

Product Datasheet

Human HCG protein (orb391476)

Description

Recombinant human HCG protein is produced

target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

CGB3; CGB5; CGB7; CGB8

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

SKEPLRPRCR PINATLAVEK EGCPVCITVN
TTICAGYCPT MTRVLQGVLP ALPQVVCNYR
DVRFESIRLP GCPRGVNPVV SYAVALSCQC
ALCRRSTTDC GGPKDHPLTC DDPRFQDSSS
SKAPPSLPS PSRLGPSDT PILPQVDHHH HHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID
P0DN86, P0DN87
Entrez
93659, 1082, 94115, 94027
Expiration Date

6 months from date of receipt.