

Product Datasheet

Human CLEC10A protein (orb391484)

Description

Recombinant human CLEC10A protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

CLEC10A

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QNSKFQRDLV TLRTDFSNT SNTVAEIQAL
TSQGSSLEET IASLKAEVEG FKQERQAGVS
ELQEHTTQKA HLGHCPCPS VCVPHSEML
LRVQQLVQDL KKLTCQVATL NNNASTEGTC
CPVNWVEHQD SCYWFSHSGM SWAEAEKYCQ
LKNAHLVIN SREEQNFVQK YLGSAYTWMG
LSDPEGAWKW VDGTDYATGF QNWKPGQPDD
WQGHGLGGGE DCAHFHPDGR
WNDDVCQRPY HWVCEAGLGQ TSQESHVDHH
HHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q8IUN9](#)

Entrez

[10462](#)