

Product Datasheet

Human CTRB1 protein (orb391564)

Description

Recombinant human CTRB1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

CTRB1

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

CGVPAIHPVL SGLSRIVNGE DAVPGSWPWQ
VSLQDKTGFH FCGGSLISED WVVTAHCGV
RTSDVVVAGE FDQGSDEENI QVLKIAKVKF
NPKFSILTVN NDITLLKLAT PARFSQTVSA
VCLPSADDDF PAGTLCATTG WGKTKYNANK
TPDKLQQAAL PLLSNAECKK SWGRRITDVM
ICAGASGVSS CMGDSGGPLV CQKDGAUTLV
GIVSWGSDTC STSSPGVYAR VTKLIPWVQK
ILAAVNDHHH HHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P17538

Entrez

1504

Expiration Date

6 months from date of receipt.