

Product Datasheet

Human DCBLD2 protein (orb391609)

Description

Recombinant human DCBLD2 protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

DCBLD2

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QQGDGCGHTV LGPESGTLTS INYPQTYPNS
TVCEWEIRVK MGERVRIKFG DFDIEDSDSC
HFNYLRIYNG IGVSRTEIGK YCGLGLQMNH
SIESKGNEIT LLFMSGIHVS GRGFLASYSV
IDKQDLITCL DTASNLEPE FSKYCPAGCL
LPFAEISGTI PHGYRDSSPL CMAGVHAGVV
SNTLGGQISV VISKGIPYYE SSLANNVTSV
VGHLSTSLFT FKTS GCYGT L GMESGVIADP
QITASSVLEW TDHTGQENSW KPKKARLKKP
GPPWAAFATD EYQWLQIDLN KEKKITGIIT
TGSTMVEHNY YVSAYRILYS DDGQKWTVYR
EPGVEQDKIF QGNKDYHQDV RNNFLPPIA
RFIRVNPTQW QOKIAMKMEL LGCQFIPKGR
PPKLTQPPPP RNSNDLKNTT APPKIAKGRA
PKFTQPLQPR SSNEFPAQTE QTTASPDIRN
TTVTPNVT KD VAVDHHHHHH

Purity

Greater than 95% as determined by

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558