

Product Datasheet

Human DFFA protein (orb391618)

Description

Recombinant human DFFA protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

DFFA

Preservatives

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MEVTGDAGVP
ESGEIRTLKP CLLRRNYSRE QHGVAASCLE
DLRSKACDIL AIDKSLTPVT LVLAEDGTIV
DDDDYFLCLP SNTKFVALAS NEKWAYNNSD
GGTAWISQES FDVDETDSGA GLKWKNVARQ
LKEDLSSIIL LSEEDLQMLV DAPCSDLAQE
LRQSCATVQR LQHTLQQVLD QREEVRQSKQ
LLQLYLQALE KEGSLLSKQE ESKAAFGEV
DAVDTGISRE TSSDVALASH ILTALREKQA
PELSLSSQDL ELVTKEPKA LAVALNWDIK
KTETVQEACE WELALRLQQT QSLHSLRSIS
ASKASPPGDL QNPKRARQDP T

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

000273

Entrez

1676