

Product Datasheet

Human DLG4 protein (orb391624)

Description

Recombinant human DLG4 protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

DLG4

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MNHKVVHHHHH HMDCLCIVTT KKYRYQDEDT
PPLEHSPAHL PNQANSPPVI VNTDTLEAPG
YELQVNGTEG EMEYEEITL RGNSGLGFSI
AGGTDNPHIG DDPSIFITKI IPGGAAAQDG
RLRVNDSILF VNEVDVREVT HSAAVEALKE
AGSIVRLYVM RRPKPAEKVM EIKLIKPKG
LGFSIAGGVG NQHIPGDNSI YVTKIIEGGA
AHKDGRLQIG DKILAVNSVG LEDVMHEDAV
AALKNTYDVV YLKVAKPSNA YLSDSYAPPD
ITTSYSQHLD NEISHSSYLG TDYPTAMTPT
SPRRYSPVAK DLLGEEDIPR EPRRIVHRG
STGLGFNIVG GEDGEGIFIS FILAGGPADL
SGELRKGDQI LSVNGVDLRN ASHEQAAIAL
KNAGQTVTII AQYKPEEYSR FEAkihDLRE
QLMNSSLGSG TASLRNPKR GFYIRALFDY
DKTKDCGFLS QALSRFGDV LHVIDAGDEE
WWQARRVHSD SETDDIGFIP SKRRVERREW
SRI KAKDWGS SSGSOGREDS VI SYETVTOM

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558